

## **MAINTENANCE AND CARE OF EQUIPMENT**

### **DO NOT STORE:**

- In **continuous heat** over 200 degrees Fahrenheit
- In **direct sunlight** for prolonged periods of time---may bleach out webbing color and webbing strength over time---(**Spectra** unaffected).
- While **damp**

## **DECONTAMINATION**

**Decontamination** is recommended after each use according to normal protocol. Any hospital-grade, E.P.A. Registered disinfectant, such as **Special Respirator Cleaner Plus** of Georgia Steel and Chemical Company, as well as **household chlorine bleach in standard dilution may be used**. Immersion time should be a minimum of **ten minutes** in either disinfecting solution. **Bleach may affect webbing color and Velcro strength in Aquatic model.**

## **CLEANING INSTRUCTIONS:**

- Use only **mild detergents** (such as Woolite) to wash
- **Cold** water wash ONLY
- **Hand Wash or Use Gentle Cycle** (preferred)
- **Use of power washer is permissible when highly soiled**
- ***HANG TO AIR-DRY ONLY***